F&P *myAIRVO 2*

User Manual

SECTION

English	A
Español	B
Français	C
繁體中文 (Chinese Traditional)	D
Melayu	E
Arabic	F

BEFORE YOU START

- This User Manual is intended for patients and healthcare professionals.
- Read this User Manual including all warnings. Failure to do so may result in injury. Keep in a safe place for future reference.
- Before the myAIRVO 2 is used for the first time, it must be set up according to the instructions in the myAIRVO 2 Technical Manual. This should be carried out by a healthcare professional or medical technician. The myAIRVO 2 needs special precautions regarding electromagnetic compliance (EMC) therefore must be installed and put into service according to the EMC information provided in this User Manual and the Technical Manual.

OTHER REFERENCES

- Refer to the myAIRVO 2 User Manual for detailed instructions for use.
- Refer to all relevant accessory User Instructions.
- Watch the training videos on the AIRVO 2 website www.fphcare.com/airvo
- For troubleshooting information, please refer to the myAIRVO 2 Technical Manual.
- Download the AIRVO 2 Simulator App to learn how to use the myAIRVO 2. You can change settings, simulate faults and test your skills. Available from the [Apple](#), [Google Play](#) and [Windows App](#) stores.
- Visit the Fisher & Paykel education & resources website at www.fphcare.com/education, to find self-paced online courses and local training events.
- If the unit is ever used by multiple patients, the unit must be cleaned and disinfected between patients according to instructions in the Disinfection Kit Manual (900PT600).
- For further assistance, please contact your Fisher & Paykel Healthcare representative.

TABLE OF CONTENTS

1. Overview	A - 2
Intended Use	A - 2
Warnings	A - 2
myAIRVO 2 and Accessories	A - 3
2. Setting up myAIRVO 2	A - 4
3. Using myAIRVO 2	A - 7
Advanced settings	A - 8
Oxygen	A - 10
Alarms	A - 11
4. Cleaning and maintenance	A - 13
Daily cleaning instructions	A - 13
Weekly cleaning instructions	A - 13
Schedule for changing accessories	A - 14
Filter replacement	A - 14
Servicing	A - 14
5. Technical Information	A - 15

1. OVERVIEW

The myAIRVO 2 is a humidifier with integrated flow generator that delivers warmed and humidified respiratory gases to spontaneously breathing patients through a variety of patient interfaces.

INTENDED USE

The myAIRVO 2 is for the treatment of spontaneously breathing patients who would benefit from receiving high flow warmed and humidified respiratory gases. This includes patients who have had upper airways bypassed. The flow may be from 2 - 60L/min depending on the patient interface. The myAIRVO 2 is for patients in homes and long-term care facilities.

USA Federal Law restricts this unit for sale by or on the order of a physician.

WARNINGS

- The unit is not intended for life support.
- Appropriate patient monitoring must be used at all times. Loss of therapy will occur if power is lost.
- Nasal delivery of respiratory gases generates flow-dependent positive airway pressure (PAP). This must be taken into account where PAP could have adverse effects on a patient.

To avoid burns:

- Use only interfaces, water chambers and breathing tubes specified in this user manual.
- Do not use accessories beyond the maximum periods of use specified in this manual.
- Before using oxygen with the unit, read all warnings in the “Oxygen” section of this manual.
- Never operate the unit if:
 - the heated breathing tube has been damaged with holes, tears or kinks,
 - it is not working properly,
 - the case screws have ever been loosened.
- Do not block the flow of the air through the unit and breathing tube.
- Locate the unit in a position where ventilation around the unit is not restricted.
- Never block the air openings of the unit or place it on a soft surface such as a bed or couch/sofa, where the filter area may be blocked. Keep the air openings free of lint, hair etc.

To avoid electric shock:

- Do not store or use the unit where it can fall or be pulled into water. If water has entered the unit enclosure, disconnect the power cord and discontinue use.
- Never operate the unit if:
 - it has been dropped or damaged,
 - it has a damaged power cord or plug,
 - it has been dropped into water.
- Avoid unnecessary removal of the power cord from the rear of the device. If removal is necessary, hold the connector during removal. Avoid pulling on the power cord.
- Return the unit to an authorized service center for examination and repair, except as outlined in this manual.

To avoid choking, or inhalation of a foreign object:

- Ensure an air filter is fitted when operating your unit.
- Never drop or insert any object into any opening or tube.


Miscellaneous:

- Prior to each patient use, ensure that the auditory alarm signal is audible by conducting the alarm system functionality check described in the Alarms section.
- Humidity output will be compromised below 18°C (64°F) and above 28°C (82°F).
- To prevent disconnection during use, especially during ambulatory use, use only heated breathing tubes specified in this manual.
- The unit is not suitable for use in the presence of a flammable, anesthetic mixture with air or oxygen or nitrous oxide.
- The myAIRVO 2 is not a sealed system. Follow hospital infection control guidelines to reduce risk of cross-contamination
- Use of accessories or power cables not specified by Fisher & Paykel Healthcare could result in increased electromagnetic emissions, decreased electromagnetic immunity and/or improper operation.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

California residents please be advised of the following, pursuant to Proposition 65:

- This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

myAIRVO 2 AND ACCESSORIES


		Nasal Cannula (20-pack)													
				Optiflow™ Junior		Optiflow™+				Optiflow™					
		Quantity per pack	Contains MR290 chamber	OPT316 (infant)	OPT318 (pediatric)	OPT942 (small)	OPT944 (medium)	OPT946 (large)	OPT970 (direct trache)	OPT980 (mask adapter)	OPT842 (small)	OPT844 (medium)	OPT846 (large)	OPT870 (direct trache)	RT013 (mask adapter)
Tube & Chamber kits	900PT500E	1				•	•	•	•	•	•	•	•	•	•
	900PT500	10				•	•	•	•	•	•	•	•	•	•
	900PT531	10	•	•	•										
	AirSpiral™														
Chamber kits	900PT560E	1		•	•	•	•	•	•	•					
	900PT560	10		•	•	•	•	•	•	•					
	900PT561	10	•	•	•	•	•	•	•	•					
	900PT290E	1	•												
	HC360	1													

Note: -E suffix denotes 1-pack

Miscellaneous	
OPT012	Wigglepads (Optiflow Junior) (20-pack)
OPT014	Oxygen Tubing (Optiflow Junior)

Some products may not be available in your country. Please contact your local Fisher and Paykel Healthcare representative.


2. SETTING UP myAIRVO 2

1. BEFORE YOU BEGIN

Place the unit on a low shelf or near the floor beside your bed. It must be placed below head height and flat. Position the device so the power cord connection to the power supply is easily accessible and able to be disconnected.

2. INSTALL WATER CHAMBER

IF USING A HC360 REUSABLE WATER CHAMBER:

With the aid of the supplied funnel, fill the chamber with enough distilled water for the period of use, but never above the 560 mL maximum water level line.

Flow setting vs usage time (HC360 reusable water chamber, at 37 °C target temperature)													
L/min	2	5	10	15	20	25	30	35	40	45	50	55	60
hrs	106	42	21	14	10	8	7	6	5	4	4	3	3


⚠️ WARNINGS

To avoid burns:

- Do not fill the water chamber with hot water.

To avoid electric shock:

- Always remove the water chamber to fill it and always fill with enough water to prevent it running out.


Fit the water chamber to the unit by pressing down the finger guard and sliding the chamber on, carefully aligning with the blue chamber port ends.

Push the chamber on firmly until the finger guard clicks into place.

Go to Step 3, "Install Heated Breathing Tube", below.

⚠️ CAUTION


Adding substances other than water can adversely affect the humidifier and delivered therapy.


IF USING THE MR290 AUTO-FILL WATER CHAMBER:

Remove the blue port caps from the chamber by pulling the tear tab upwards then remove the bracket holding the water supply tube. Fit the supplied adapter over the two vertical ports on the chamber and push on fully then clip the water supply tube into position.

Fit the MR290 chamber as described above for the HC360 chamber.

Attach the water bag to the hanging bracket at least 20 cm (8") above the unit, and push the bag spike into the fitting at the bottom of the bag. Open the vent cap on the side of the bag spike. The chamber will now automatically fill to the required level and maintain that level until the water bag is empty. Use only distilled water and always ensure sufficient water is in the water bag to prevent it from running out.

Check that water flows into the chamber and is maintained below the maximum water level line. If the water level rises above the maximum water level line, replace the chamber immediately.

**MR290: Flow setting vs usage time
(Water bag 900PT401, 1000 mL, at 37 °C target temperature)**

L/min	2	5	10	15	20	25	30	35	40	45	50	55	60
hrs	189	75	37	25	18	15	12	10	9	8	7	6	6

⚠️ WARNINGS

To avoid burns:

- Do not start the unit without the water chamber in place.
- Do not touch the heater plate, water chamber or chamber base during use.
- The water in the chamber becomes hot during use. Exercise caution when removing and emptying the chamber.


To avoid electric shock:

- When handling the unit with the water chamber in place, avoid tilting the machine to prevent any chance of water entering the unit enclosure.
- Empty all the water from the water chamber before transporting the unit.

⚠️ CAUTIONS

To ensure optimal therapy (MR290 only):

- Do not use the auto-fill MR290 chamber if it has been dropped or been allowed to run dry, this could lead to the chamber over filling.
- Do not use the MR290 chamber if the water level rises above the maximum water level line as this may lead to water entering the patient's airway.

3. INSTALL HEATED BREATHING TUBE

One end of the heated breathing tube has a blue plastic sleeve. Lift the sleeve and slide the connector onto the unit. Push the sleeve down to lock.

⚠️ WARNINGS

To avoid burns:

- Do not modify the breathing tube or interface in any way.
- Do not allow the breathing tube to remain in direct contact with skin for prolonged periods of time. The healthcare professional shall assess the conditions for safe contact, such as duration and skin condition.
- Do not add heat above ambient levels to any part of the breathing tube or interface e.g. by covering with a blanket or by heating with infrared radiation, an overhead heater, or an incubator.
- Do not use an insulating sleeve or any similar accessories which are not recommended by Fisher & Paykel Healthcare.

⚠️ CAUTIONS

- Position the heated breathing tube away from any electrical monitoring leads (EEG, ECG/EKG, EMG, etc), to minimize any possible interference with the monitored signal.

CONDENSATE MANAGEMENT

The unit must be placed below head height and flat, this allows condensate to drain towards the water chamber, away from the patient.


If excess condensate accumulates in the heated breathing tube, disconnect the patient interface from the heated breathing tube, drain the condensate by lifting the patient end of the tube, allowing the condensate to run into the water chamber.

At higher target flow rates, it may be necessary to first reduce the target flow rate to 30 L/min or below, to ensure the condensate drains into the water chamber.

Minimize local sources of cooling acting on the heated breathing tube, such as a fan to cool the patient or an air-conditioning unit/vent

If condensate persists, consider turning the target temperature down. A lower target temperature will decrease the humidity output of the unit, decreasing the level of condensation.

Note: The temperature and humidity level delivered to the patient will also be reduced.




4. SELECT PATIENT INTERFACE

The myAIRVO 2 can be used with a variety of patient interfaces. Read the separate user instructions for the patient interface that will be used, including all warnings.

⚠️ WARNINGS


To avoid burns:

- Do not modify the breathing tube or interface in any way.
- Do not use any patient interfaces not listed here.

Nasal cannula			Tracheostomy interface	Mask interface adapter	
					
Optiflow™+ OPT942E OPT944E OPT946E	Optiflow™ OPT842E OPT844E OPT846E	Optiflow™ Junior/ OPT316/OPT318 (Refer to "Using AIRVO 2" - "Junior Mode")	OPT970E / OPT870E	OPT980E / RT013E (with mask) Note that the OPT980E/RT013E Mask Interface Adapter is designed to be used with vented masks only. Do not use sealed masks.	

All patient interfaces are Type BF applied parts.


The following table shows the target dew-point temperature settings and target flow settings able to be used with these interfaces.

Patient Interface	°C			L/min											
	31	34	37	2	5	10	15	20	25	50	55	60	
 OPT316	●	●	●	2		20									
 OPT318	●	●	●	2		25									
 OPT942E (S)	●	●	●			10		50							
 OPT944E (M)	●	●	●			10		60							
 OPT946E (L)	●	●	●			10		60							
 OPT970E	●	●	●			10		60							
 OPT980E	●	●	●			10		60							
 OPT842E (S)	●	●	●			10		50							
 OPT844E (M)	●	●	●			10		60							
 OPT846E (L)	●	●	●			10		60							
 OPT870E	●	●	●			10		60							
 RT013E	●	●	●			10		60							

Low temperature ambient conditions may prevent the unit from reaching a 37 °C target temperature setting at high target flow settings. In these cases, consider decreasing the target flow setting.

At altitude, the maximum flow rates achievable may be lower than those in the above table, by approximately 5 L/min per 1000 m (3000 ft).

3. USING myAIRVO 2

1. SWITCH ON UNIT

Plug the unit's power cord into the mains/utility power socket. The connector at the other end of the power cord should be well secured to the rear of the unit.


⚠️ WARNINGS

To avoid electric shock:

- Ensure that the unit is dry before plugging into the mains/utility power socket.

Switch on the unit by pressing the On/Off button for 5 seconds.


2. WARM-UP

The unit will begin to warm up. You will see a warm-up symbol on the screen.

"Warm-up" symbol


3. CONNECT THE PATIENT INTERFACE

"Ready for use" symbol

When the "Ready for Use" symbol appears on the display, connect the patient interface to the heated breathing tube. When you first use the unit, the air will feel warm. Continue to breathe normally.


4. AFTER USE

When you have finished using the unit, remove your interface and drain any excess condensate in the breathing tube by lifting the patient end of the tube, and allowing the condensate to run into the water chamber.

5. DRYING MODE

Then press and hold the On/Off button for 3 seconds until a melody sounds. The unit will automatically enter Drying Mode and dry the tube so it is ready for you to use next time. Drying Mode runs for 99 minutes. The unit will automatically turn off when it is finished.

⚠️ WARNINGS

To avoid burns:

- Do not wear the interface during Drying Mode. The air is hot and dry and may cause injury.
- Do not remove the water chamber until drying mode has been completed.


To switch the unit off without completing Drying Mode (this is not recommended), hold down the On/Off button for 5 seconds.

If you unplug the unit's power cord from the mains/utility power socket while the unit is still running, the "Power Out" alarm will sound. Press the "Audio Pause" button to silence this alarm.

ADVANCED SETTINGS

When you see the “Warm-up” or “Ready for use” symbols, you can press the Mode button to view and change advanced settings.

TARGET DEW-POINT TEMPERATURE

You can set the myAIRVO 2 to three target dew-point temperature settings:

- 37°C (98.6°F)
- 34°C (93°F) [if compliance at 37°C is a problem]
- 31°C (88°F) [for face masks only].

You may not have access to all settings, if:

- the unit is in Junior Mode (limited to 34 °C),
- the unit was initially set up with tighter limits.


The myAIRVO 2 will remember its target dew-point temperature setting when you switch it off.

To change the target dew-point temperature setting:

Press the Up and Down buttons to choose the new setting.

The large number in the center of the screen shows your chosen setting.

The small numbers near the arrow show the minimum and maximum accessible settings.

Press the Mode button to move on to the next screen.

TARGET FLOW

You can set the myAIRVO 2 to flows between 10 L/min and 60 L/min, in increments of 1 L/min (10-25 L/min) and 5 L/min (25-60 L/min).

You may not have access to all settings, if:

- the unit is in Junior Mode (limited to 2 - 25 L/min, in increments of 1 L/min)
- the unit was initially set up with tighter limits.

The myAIRVO 2 will remember its target flow setting when you switch it off.

At altitude, the maximum flow rates achievable may be lower than those listed above, by 5 L/min per 1000 m (3000 ft).

To change the target flow setting:


Press the Up and Down buttons to choose the new setting.

The large number in the center of the screen shows your chosen setting.

The small numbers near the arrow show the minimum and maximum accessible settings.

Press the Mode button to move on to the next screen.


DAY/NIGHT MODES

You can set the myAIRVO 2 to “Day” mode or “Night” mode.

In “Night” mode, some of the myAIRVO 2 sounds will be made quieter. The display will become dimmed. Alarms will be unaffected.

The myAIRVO 2 will remember its Day/Night setting when you switch it off.

To change the Day/Night setting:

Press the Up and Down buttons to choose the new setting.


→ “Day”

→ “Night”

Press the Mode button to move on to the next screen.

COMPLIANCE

This screen displays three pieces of compliance data:

<i>Total hours used</i>	Displays the total number of hours that the unit has been switched on.
<i>Hours per day</i>	Displays the average number of hours that the unit has been used per day.
<i>Checksum</i>	Displays usage information for the medical practitioner.

Press the Mode button to return to the “Warm-up”/“Ready for use” screen.

JUNIOR MODE

If the patient will be using an Optiflow Junior nasal cannula (OPT316/OPT318), you must activate Junior Mode. Do not use Junior mode for other patient interfaces.


Junior Mode limits the target settings to: 34 °C and 2 - 25 L/min, in increments of 1 L/min.

To activate Junior Mode:


You must be able to see the “Warm-up” symbol or the “Ready for use” symbol to activate Junior Mode.

Hold the Mode button for 5 seconds.

New target settings


The target settings for dew-point temperature and flow will be changed automatically. The colorful icons in the corners of the screen indicate that this unit is in Junior Mode.

To deactivate Junior Mode, follow the same procedure: hold the Mode button for 5 seconds.

If you are unable to activate Junior Mode, it is possible that Junior Mode may not have been enabled for your device. Contact your Fisher & Paykel Healthcare representative.

OXYGEN

You can connect up to 15 L/min of supplementary oxygen from a regulated supply to the myAIRVO 2. Connect the output from the oxygen source to the oxygen inlet port on the side of the unit. Make sure you push the oxygen tube firmly onto this connection port.

The fraction of oxygen you breathe with this air/oxygen mixture is determined by the airflow setting on the unit and the oxygen flow connected to the unit's oxygen inlet port.

The following table gives the approximate oxygen fraction delivered for the range of unit and oxygen airflows (at sea level). The oxygen fractions given assume that the oxygen source is a home oxygen concentrator. These values will be higher if the oxygen source is bottled oxygen. At flows less than 10 L/min, the oxygen fraction delivered varies significantly with small changes in input oxygen flow. Oxygen flow settings should be titrated according to blood saturation levels.

		myAIRVO 2 Target Flow Setting (L/min)										
		10	15	20	25	30	35	40	45	50	55	60
Oxygen Flow (L/min)	1	29	27	25	24	24	23	23	23	23	23	22
	2	38	32	29	28	26	26	25	25	24	24	24
	3	45	37	33	31	29	28	27	26	26	25	25
	4	53	42	37	34	32	30	29	28	27	27	26
	5	60	48	41	37	34	33	31	30	29	29	28
	7	75	58	50	44	40	37	35	34	32	31	31
10	93	74	61	54	49	45	42	39	37	36	35	

It is important that the physician prescribing your oxygen therapy approves both the flow and oxygen settings and that you do not adjust these prescribed settings without consulting them.

Check that suitable blood saturation levels are achieved at the prescribed flow.

Use continuous oxygen monitoring on patients who would desaturate significantly in the event of disruption to their oxygen supply.

⚠️ WARNINGS

Before using oxygen with the unit, read all of the following warnings:


- The use of oxygen requires that special care be taken to reduce the risk of fire. Accordingly, for safety it is necessary that all sources of ignition (e.g. electrocautery or electrosurgery) be kept away from the unit and preferably out of the room in which it is being used. Oxygen should not be used while smoking or in the presence of an open flame. The unit should be located in a position where ventilation around the unit is not restricted.
- A spontaneous and violent ignition may occur if oil, grease or greasy substances come in contact with oxygen under pressure. These substances must be kept away from all oxygen equipment.
- Ensure that the myAIRVO 2 is switched on before connecting oxygen.
- Oxygen must only be added through the special oxygen inlet port on the back of the unit. To ensure that oxygen enters the unit correctly, the oxygen inlet port must be fitted properly to the filter holder and the filter holder must be fitted properly to the unit. The power cord connector should also be well secured.
- Ensure that the myAIRVO 2 target flow rate is higher than the supplementary oxygen flow rate, to avoid excess oxygen being vented into the surroundings.
- Do not connect more than 15 L/min of supplementary oxygen to the myAIRVO 2.
- The oxygen concentration delivered to the patient can be affected by changes to the flow setting, oxygen setting, patient interface or if the airpath is obstructed.

When finished, turn off the oxygen source. Remove the output of the oxygen source from the oxygen inlet port on the back of the unit.

⚠️ WARNINGS

To avoid burns:

- The oxygen flow must be turned off when the unit is not operating, so that oxygen does not build up inside the device.

ALARMS

The myAIRVO 2 has visual and auditory alarms to warn you about interruptions to your treatment. These alarms are generated by an intelligent alarm system, which processes information from the sensors and target settings of the unit and compares this information to pre-programmed limits.

ALARM SIGNALS

	Symbols	Meaning
Visual alarm signal		
		Alarm condition.
		Audio paused.
Auditory alarm signal		
3 beeps in 3 seconds. Repeated every 5 seconds.		Press this button to mute the auditory alarm for 115 seconds. The auditory alarm can be reactivated by pressing this button again.

ALARM CONDITIONS

All of the alarms listed below have been assessed as “Medium Priority”. These priorities have been allocated for an operator’s position within 1 meter of the device. The unit also uses an internal priority-ranking system. If multiple alarm conditions occur simultaneously, the unit will display the highest-priority alarm.

The following table lists all of the alarm conditions from highest-priority to lowest priority, their causes, possible solutions and delays. Alarm conditions that affect oxygen delivery require an immediate response to assess the patient’s saturation levels. Alarm conditions that affect humidity delivery require a prompt response to assess potential drying of mucus and associated blockages.

The following alarm delays assume operation in ‘Ready for use’ mode.

Message	Meaning	Affects delivery of:	Delays
<i>Fault (E###)</i>	<i>The unit has detected an internal fault and has shut itself down.</i> Switch the unit off and then restart. If the problem persists, note the fault code and contact your Fisher & Paykel Healthcare representative.	Oxygen, humidity.	< 5 seconds
<i>Check tube</i>	<i>The unit cannot detect the heated breathing tube.</i> Check that the heated breathing tube is not damaged and that it is plugged in correctly. If the problem persists, then change the heated breathing tube.	Oxygen, humidity.	< 5 seconds
<i>Check for leaks</i>	<i>The unit has detected a leak in the system.</i> The most likely cause is that the water chamber has been removed or has not been pushed into place correctly. Check that the heated breathing tube is not damaged and that it is plugged in correctly. Check that the nasal interface is fitted. Check that the filter is fitted.	Oxygen, humidity.	< 120 seconds
<i>Check for blockages</i>	<i>The unit has detected a blockage in the system.</i> Check the heated breathing tube or patient interface for blockage. Check the air filter and filter holder for blockage. Check whether the unit should be in Junior Mode. If the patient will be using an Optiflow Junior nasal cannula (OPT316/OPT318), you must activate Junior Mode.	Oxygen, humidity.	< 10 seconds
<i>O₂ too low</i>	<i>The measured oxygen level has fallen below the allowed limit.</i> Check that the oxygen source is still operational and is correctly connected. Adjust the level of oxygen from the oxygen source as necessary.	Oxygen	< 20 seconds
<i>O₂ too high</i>	<i>The measured oxygen level has exceeded the allowed limit.</i> Check that the myAIRVO flow rate is set correctly. Adjust the level of oxygen from the oxygen source as necessary.	Oxygen	< 20 seconds
<i>Check operating conditions</i>	<i>The unit has detected that it is operating in unsuitable ambient conditions.</i> This alarm may be caused by a sudden change in ambient conditions. Leave the unit running for 30 minutes. Switch the unit off and then restart	Humidity	60 +/- 6 seconds

(continued)			
Message	Meaning	Affects delivery of:	Delays
Cannot reach target flow	<p>The unit cannot reach the target flow setting.</p> <p>Check the heated breathing tube or patient interface for blockage. Check whether the target flow setting is too high for the patient interface being used (refer to "Setting up myAIRVO 2" - "Select Patient Interface"). You will be prompted for acknowledgement.</p> <p>⚠️ WARNINGS</p> <ul style="list-style-type: none"> The oxygen concentration delivered to the patient can be affected by changes to the flow setting. Adjust the level of oxygen from the oxygen source as necessary. 	Oxygen	< 120 seconds
Check water	<p>The chamber has run out of water.</p> <p>If using the HC360 reusable chamber: Remove the chamber and refill.</p> <p>If using the MR290 auto-fill chamber: When a chamber runs dry, the chamber float may be damaged. Replace the chamber and water bag.</p> <p>To ensure continual humidification, always ensure that the water chamber and/or water bag are not allowed to run out of water.</p>	Humidity	< 30 minutes
Cannot reach target temperature	<p>The unit cannot reach the target temperature setting.</p> <p>You will be prompted for acknowledgement.</p> <p>The most likely cause for this is that the unit is operating at a high flow rate in low ambient conditions. Consider increasing ambient conditions to match the recommended operating conditions or decreasing the target flow setting.</p> <p>⚠️ WARNINGS</p> <ul style="list-style-type: none"> The oxygen concentration delivered to the patient can be affected by changes to the flow setting. Adjust the level of oxygen from the oxygen source as necessary. 	Humidity	30 +/- 3 minutes
[Power out]	<p>The unit has been disconnected from the mains/utility power socket.</p> <p>No visual alarm. The auditory alarm will sound for 120 seconds. If power is reconnected in this time, the unit will automatically restart.*</p> <p>⚠️ WARNINGS</p> <ul style="list-style-type: none"> Appropriate patient monitoring must be used at all times. Loss of therapy will occur if power is lost. 	Oxygen, humidity.	< 5 seconds

ALARM LIMITS

Most alarm limits are pre-programmed. The exceptions are listed below. These alarm limits may be changed to other values by authorized personnel. Changes will be preserved during or after any power loss.

Alarm condition	Factory-set alarm limit	Possible preset values
O ₂ too low	21% O ₂	21 or 25% O ₂
O ₂ too high	90% O ₂	30 - 95% O ₂ in 5% increments

⚠️ WARNINGS

- A hazard can exist if different alarm presets are used on different units within any single area, eg. long term care facility
- Alarm limits set to extreme values can render the alarm system useless.

CHECKING ALARM SYSTEM FUNCTIONALITY

The functionality of the alarm system can be checked at any time when the unit is turned on. Remove the heated breathing tube. You should see the "Check tube" visual alarm signal and hear the auditory alarm signal. If either alarm signal is absent, do not use the unit. Contact your Fisher & Paykel Healthcare representative.

AUDITORY INFORMATION SIGNALS

In addition to auditory alarm signals, auditory information signals are provided. These are described below.

Melody	Meaning
Ascending sequence of 5 tones	The "Ready for use" symbol has appeared
Ascending sequence of 3 tones	Activation/deactivation of Junior Mode
Decending scale of 3 tones (within 2 seconds)	Drying Mode has been activated
Single tone every 5 seconds	Measured oxygen level \geq 33% at turn-off
Single tone every 30 seconds	Measured oxygen level $>$ 95%

4. CLEANING AND MAINTENANCE

It is important to carefully follow the instructions in this section, to keep the device clean and safe for use and to extend the life of the consumables.

The following instructions are for single-patient home use. If the unit is ever used by multiple patients, the unit must be cleaned and disinfected between patients according to instructions in the Disinfection Kit Manual (900PT600). Single-patient use accessories must be disposed of between patients to prevent cross-contamination.

Standard aseptic techniques to minimize contamination should be followed when handling the unit and accessories. This includes proper hand-washing, avoiding hand contact with connection ports, safe disposal of the used consumables and suitable storage of the unit after cleaning and disinfection.

DAILY CLEANING INSTRUCTIONS

Run Drying Mode / Rinse the patient interface and water chamber

1. Allow Drying Mode to run after use (refer “Using myAIRVO 2” - “Drying Mode”).
2. Remove the interface, rinse in drinking-quality water then reconnect to the heated breathing tube whilst still in Drying Mode to dry the interface.
- 3.

If using the HC360 reusable chamber:

After Drying Mode is complete, remove the water chamber by pushing down the finger guard and pulling out the chamber. Wash and rinse the chamber then refill it with sufficient distilled water for the next use.

If using the MR290 chamber:

Do not wash or remove this chamber.

WEEKLY CLEANING INSTRUCTIONS

Clean the patient interface, water chamber and myAIRVO 2

1. Switch off the unit and unplug from the mains/utility power socket.
2. Remove the heated breathing tube and drain any excess condensate.
3. Remove the interface from the heated breathing tube, wash it in warm water with mild dishwashing detergent added, rinse it in drinking-quality water, then reconnect it to the heated breathing tube.
4. Remove the water chamber.
5. **If using the HC360 chamber:** Pour out and discard the remaining water. Remove the chamber base. Wash the chamber top and base in mild dishwashing detergent then rinse. Soak the chamber in a solution of vinegar (1 part) and water (2 parts) for 10 minutes. Rinse and dry.
If using the MR290 chamber: Do not wash this chamber. Carefully put the MR290 chamber aside.
6. Thoroughly wipe the inside of the heated breathing tube connection port with a clean, low-lint cloth dipped in warm water with mild dishwashing detergent added.
7. Wipe the exterior of the unit with a clean, damp (not wet) cloth dipped in warm water with mild dishwashing detergent added. Do not use harsh abrasives or solvents, as these may damage the unit.
8. Refit the heated breathing tube.
9. **If using the HC360 chamber:** Refit the chamber.
If using the MR290 chamber: Refit the MR290 chamber and reconnect to the water bag. Check that water flows into the chamber and is maintained below the maximum water level line. If the water level rises above the maximum water level line, replace the chamber immediately.
10. Reconnect the unit to the mains/utility power socket.
11. The unit is now ready for another week of use.

SCHEDULE FOR CHANGING ACCESSORIES



The accessories for the unit must be changed frequently to avoid the risk of infection. Parts should be replaced immediately if they are damaged or discolored; otherwise they must be replaced within the periods shown in the following table. These periods assume that the correct daily and weekly cleaning procedures and maintenance schedule described above are adhered to. If these procedures and schedules are not followed, the maximum periods of use will change to those stated in the AIRVO 2 manual. These accessories are for single-patient use only.

Maximum period of use	Part number and description
1 week	Optiflow Junior interfaces OPT316 Nasal Cannula - Infant OPT318 Nasal Cannula - Pediatric
1 month	All other patient interfaces OPT942 / OPT942E Optiflow™+ Nasal Cannula - Small OPT944 / OPT944E Optiflow™+ Nasal Cannula - Medium OPT946 / OPT946E Optiflow™+ Nasal Cannula - Large OPT970 / OPT970E Tracheostomy Interface OPT980 / OPT980E Mask Interface Adapter OPT842 / OPT842E Optiflow™ Nasal Cannula - Small OPT844 / OPT844E Optiflow™ Nasal Cannula - Medium OPT846 / OPT846E Optiflow™ Nasal Cannula - Large OPT870 / OPT870E Tracheostomy Interface RT013 / RT013E Mask Interface Adapter - 22mm
2 months	All tube & chamber kits 900PT561 AirSpiral™ heated breathing tube, MR290 auto-fill chamber and adapter 900PT560 / 900PT560E AirSpiral™ heated breathing tube 900PT290E MR290 auto-fill chamber and adapter 900PT500 / 900PT500E Heated breathing tube 900PT501 Heated breathing tube, MR290 auto-fill chamber and adapter 900PT531 Junior heated breathing tube, MR290 auto-fill chamber and adapter (for use with OPT316/OPT318 only)
3 months or 1000 hours	900PT913 Air filter (or more often if significantly discolored)
Reusable	HC360 Reusable water chamber

Some products may not be available in your country. Please contact your local Fisher and Paykel Healthcare representative.

FILTER REPLACEMENT

After the myAIRVO 2 has been switched on for a total of 1000 hours, a prompt will appear indicating that an air filter change is due. Follow the steps below if filter change is due:

1. Take the filter holder from the back of the unit and remove the filter.
2. Replace the old filter with a new one.
3. Reattach the filter holder to the unit (clip the bottom of the filter holder in first, then rotate it upwards until the top clips into place).
4. Press the Mode button to move on to the “Replace now” screen.
5. Press the Up button to select “Now”.
6. Press the Mode button to confirm.
The hours counter will be reset to zero.

If you choose the “Later” option, the prompt will continue to appear whenever the unit is switched on.

SERVICING

This device contains no internal serviceable parts.

Refer to the myAIRVO 2 Technical Manual for a list of external spare parts.

5. TECHNICAL INFORMATION

SYMBOL DEFINITIONS

	For safety reasons, refer to the instructions for use		Class II equipment
	Caution		Catalogue number
	Consult instructions for use		Serial number
	Warning, hot surface		Batch code
	Manufacturer		Humidity range
	Date of manufacture		Temperature range
	Date of shelf life expiry	IP22	Protected against ingress of small objects and water drops
	Type BF applied part		EU representative
Rx only	(USA) Federal Law restricts this device to sale by, or on the order of a physician.		CE Mark
	Alarm symbol		Power on/off (standby)
	Alarm pause		Regulatory Compliance Mark (RCM)

PRODUCT SPECIFICATIONS

<i>Dimensions</i>	295 mm x 170 mm x 175 mm (11.6" x 6.7" x 6.9")	<i>Target temperature settings</i>	37, 34, 31 °C
<i>Weight</i>	2.2 kg (4.8 lb) unit only, 3.4 kg (7.5 lb) packaged in bag incl. accessories	<i>Humidity performance</i>	>33 mg/L at 37 °C target >12 mg/L at 34 °C target >12 mg/L at 31 °C target
<i>Supply frequency</i>	50-60 Hz	<i>Maximum temperature of delivered gas</i>	43 °C (109 °F) (in accordance with ISO 80601-2-74)
<i>Supply voltage/current</i>	100-115 V 2.2 A (2.4 A max†) 220-240 V 1.8 A (2.0 A max†)	<i>Maximum surface temperature of applied parts</i>	44 °C (111 °F) (in accordance with ISO 80601-2-74)
<i>Sound pressure level</i>	Alarms exceed 45dbA @ 1 m	<i>Flow range (default)</i>	10-60 L/min*
<i>Auditory alarm pause</i>	115 seconds	<i>Flow range (Junior Mode)</i>	2-25 L/min*
<i>Expected service life</i>	5 years	<i>Maximum oxygen input</i>	60 L/min
<i>Serial port</i>	The serial port is used for downloading product data, using F&P Infosmart™ software.	<i>Oxygen analyzer accuracy</i>	< ± 4 % (within the range 25-95% O ₂) Operating conditions: 18-28 °C (64-82 °F), 30-70% RH
<i>Warm-up time</i>	10 minutes to 31 °C (88 °F), 30 minutes to 37 °C (98.6 °F) using a MR290 chamber with flow rate of 35 L/min and starting temperature 23 ± 2 °C (73 ± 3 °F)		

* Flow rates are measured in BTPS (Body Temperature/Pressure, Saturated)

† Inrush current may reach 50A

OPERATING CONDITIONS

<i>Ambient temperature</i>	18 - 28 °C (64 - 82 °F)
<i>Humidity</i>	10 - 95% RH
<i>Altitude</i>	0 - 2000 m (6000 ft)
<i>Mode of operation</i>	Continuous operation

STORAGE AND TRANSPORT CONDITIONS

myAIRVO

<i>Ambient temperature</i>	-10 - 60 °C (14 - 140 °F)
<i>Humidity</i>	10 - 95% RH, non-condensing

Tube & chamber kits

<i>Ambient temperature</i>	-10 - 50 °C (14 - 122 °F)
<i>Humidity</i>	10 - 95% RH, non-condensing

The unit may require up to 24 hours to warm up or cool down from the minimum or maximum storage temperature before it is ready for use.

WARNING

- Do not use the unit at an altitude above 2000 m (6000 ft) or outside a temperature range of 18 - 28 °C (64 - 82 °F). Doing so may affect the quality of the therapy or injure the patient.

Designed to conform to the requirements of:

IEC 60601-1:2005 + A1:2012
IEC 60601-1-2:2014
ANSI/AAMI 60601-1:2005/(R) 2012
CAN/CSA-C22.2 No. 60601-1:2014
EN 60601-1:2006 + A1:2013
ISO 80601-2-74:2017

The unit complies with the electromagnetic compatibility requirements of IEC 60601-1-2. In certain circumstances, the unit may affect or be affected by nearby equipment due to the effects of electromagnetic interference. Excessive electromagnetic interference may affect the therapy delivered by the unit. If this should happen, try moving the unit or the location of the unit causing interference, or alternatively consult your healthcare provider. To avoid potential interference, do not place any part of the device or accessories within 30 cm (12") of any portable or mobile radio frequency communication equipment.

Accessory equipment connected to the serial port of the device must be certified to either IEC 60601-1 or IEC 60950-1. Furthermore all configurations shall comply with the system standard IEC 60601-1-1. Anyone who connects additional equipment to the signal input part or signal output part configures a medical system and is therefore responsible for ensuring that the system complies with the requirements of the system standard IEC 60601-1-1. If in doubt, consult the technical services department or your local representative.

DISPOSAL INSTRUCTIONS

Unit Disposal Instructions

This unit contains electronics. Please do not discard with regular waste. Return to Fisher & Paykel Healthcare or dispose according to local guidelines for disposing of electronics. Dispose according to Waste Electrical and Electronic Equipment (WEEE) directive in European Union.

Consumables Disposal Instructions

Place the interface, breathing tube and chamber in a bag at the end of use and discard with regular waste.

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For more information please contact your local Fisher & Paykel Healthcare representative

Manufacturer

Fisher & Paykel Healthcare Ltd,
15 Maurice Paykel Place,
East Tamaki, Auckland 2013

PO Box 14 348 Panmure,
Auckland 1741,
New Zealand

Tel: +64 9 574 0100
Fax: +64 9 574 0158
Email: info@fphcare.co.nz
Web: www.fphcare.com

Australia (Sponsor)

Fisher & Paykel Healthcare Pty
Ltd,
19-31 King Street, Nunawading,
Melbourne, Victoria 3131.

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Tel: +90 312 354 34 12
Fax: +90 312 354 31 01

UK

Fisher & Paykel Healthcare Ltd,
Unit 16, Cordwallis Park,
Clivemont Road, Maidenhead,
Berkshire SL6 7BU, UK

Tel: 0800 132 189
Fax: +44 1628 626 146

USA/Canada

Tel: 1800 446 3908
or +1 949 453 4000
Fax: +1 949 453 4001

F&P, AIRVO, Optiflow and AirSpiral are all trademarks of Fisher & Paykel Healthcare Ltd.
For patent information, see www.fphcare.com/ip

 0123

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Ürünleri Ticaret Limited Şirketi,